

TARDOC 2026 - Tariff positions

Reference for wound care

This reference supports physicians and wound-care centres in billing services under the new Swiss outpatient tariff.



TARDOC - PHYSICIAN SERVICES

CODE	DESCRIPTION	TP
AA.00.0060	Consultation before or after an inpatient or outpatient treatm...	3.84
MK.15.0030	Wound treatment, first dressing	21.24
MK.15.0040	Set-up and clean-up for wound treatment, per minute	4.25
MK.15.0050	Wound treatment, each additional cm	8.50
MK.15.0060	Debridement of an open wound, per 4 cm²	50.99
MK.15.0070	Set-up and clean-up for debridement, per minute	4.25
MK.15.0080	Debridement of an open wound, each additional 4 cm²	17.00



OUTPATIENT FLAT RATES

CODE	DESCRIPTION	TP
C09.12A	Skin harvesting or coverage with a full-thickness skin graft, in...	5 401.23
C09.12B	Surgical procedures on the skin, subcutaneous tissue or asso...	3 044.33
C09.12C	Wound treatment involving complex structures with imaging...	1 973.01
C09.12D	Procedures on the extremities, complete wound closure, dee...	1 206.57
C09.12E	Wound treatment and debridement with conventional radiogr...	842.87
C09.12F	Wound treatment and debridement involving complex struct...	656.05

YOUR CONTACT

+41 44 537 83 60
mail@micromedical.ch

MICRO-MEDICAL AG

Langäulistrasse 8
CH-9470 Buchs
Switzerland



TARDOC - NON-PHYSICIAN SERVICES

CODE	DESCRIPTION	TP
AK.00.0070	Non-physician treatment of a wound without complex structures, per minute -...	1.73
AK.00.0080	Non-physician wound cleaning or dressing change, per minute · (up to 10 times per...	1.73
AK.00.0090	Non-physician removal of dressings and suture material, per session · (up to 1 tim...	34.56



TARIFF TYPE

Consumables

MATERIAL WITH GTIN

Endoform® and Myriad® are medical devices with a GTIN. In outpatient TARDOC services, medical devices with a GTIN may generally be billed via tariff type 402, provided the tariff requirements are met (in particular publication of the GTIN and the applicable billing rules).



IMPORTANT NOTES

This overview is provided for guidance only and does not replace a tariff-legal review. Billing is governed by TARDOC, the outpatient flat rates and the tariff agreements in force.

Further information and application guides: micromedical.info

TARDOC AT A GLANCE



Outpatient physician services



Wound treatment, debridement, dressing changes, follow-ups



Materials with a GTIN can – provided the requirements are met – be billed separately via tariff type 402.



INDICATIONS & PRODUCTS

Endoform®

Ideal for superficial to deep chronic wounds

- ✓ Chronic wounds (diabetic foot ulcers, venous...
- ✓ Pressure ulcers
- ✓ Postoperative wounds
- ✓ Wound dehiscence
- ✓ Traumatic wounds



10 × 12,7



5 × 5



Ø 2,5

Myriad Matrix™

Ideal for deep, complex or surgical defects

- ✓ Deep soft-tissue defects
- ✓ Exposed anatomical structures
- ✓ Complex surgical defects
- ✓ Traumatic defects
- ✓ Oncological surgery



10 × 20



10 × 10



5 × 5



EVIDENCE-BASED

Extensive clinical experience in regener...



SAFE & BIOCOMP...

Decellularised extracellular matrix -...



EFFICIENT

Supports wound healing and shortens treatme...



ECONOMICAL

Transparent billing via TARDOC and tariff typ...